Glutathion HPLC Kit

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1. INTENDED USE	17
2. SUMMARY AND EXPLANATION OF THE TEST	17
3. PRINCIPLE OF THE TEST	17
4. MATERIAL SUPPLIED	19
5. MATERIAL REQUIRED BUT NOT SUPPLIED	19
6. PREPARATION AND STORAGE OF REAGENTS	20
7. PRECAUTIONS	20
8. SPECIMEN COLLECTION AND PREPARATION	20
9. ASSAY PROCEDURE	21
PROCEDURAL NOTES	21
SAMPLE AND STANDARD PREPARATION	22
CHROMATOGRAPHIC CONDITIONS	23
10. TREATMENT OF THE COLUMN	23
11. RESULTS	23
CALCULATION	23
TYPICAL CHROMATOGRAM	24
12. LIMITATIONS	25
13. QUALITY CONTROL	25
EXPECTED VALUES	25
CONTROLS	25
14. PERFORMANCE CHARACTERISTICS	26
PRECISION AND REPRODUCIBILITY	26
LINEARITY	26
DETECTION LIMIT	26
15. DISPOSAL	26
16. TROUBLESHOOTING	27
17. REFERENCES	28
18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	28

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
Glutathione HPLC Kit

For the determination of Glutathione in EDTA-blood

Gültig ab/valid from 30.06.2008

KC 1800

CAL
INT STD
CTRL 1
CTRL 2

Immundiagnostik AG, Stubenwald-Allee 8a, D 64625 Bensheim
Phone: ++49 6251 70190-0
Fax: ++ 49 6251 849430
e.mail: Info@immundiagnostik.com
www.immundiagnostik.com

1. INTENDED USE

The Immundiagnostik Assay is intended for the quantitative determination of glutathione in EDTA-blood. This assay is designed for in vitro diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Glutathione (GSH) is an intracellular tripeptide common in all tissues, which protects the cells against oxidative processes. It has important functions in several metabolic pathways like activation or inhibition of enzymes, and transport of molecules and at the transport of amino acids.

A very important function is stabilizing SH-groups in proteins and other molecules to maintain a reducing intracellular environment.

Alterations in the glutathione status are involved in the pathogenesis of several diseases:

Reperfusion damage, liver injury, cancer, diabetes mellitus, cataract, inflammatory diseases, chronic lymphatic oedema, radiation damages.

Altered glutathione concentrations might also be due to pollution, cigarette smoke, side effects of drugs and aging.

Most of the cellular glutathione is reduced (in EDTA-blood approx. 90 %), only a minor amount of 10 % is oxidized (GSSG). This steady state is maintained by the NADPH-dependent glutathione reductase.

In case of oxidative stress GSH is needed for several reactions of the primary and secondary anti oxidative protection.

3. PRINCIPLE OF THE TEST

The determination of glutathione starts by adding a dilution solution to the sample and dividing it in two aliquots.

The measurement of the reduced fraction is performed by the addition of reaction buffer and derivatisation solution. After an incubation of 20 minutes, in which GSH is transformed into a fluorescent product, a precipitation solution is added to separate higher molecular substances.

The measurement of the total glutathione is performed by the addition of the reduction solution, internal standard and derivatisation solution. After that, the sample is handled like the reduced fraction. Depending on the reduction solution all the glutathione (GSH and GSSG) is reduced to GSH.

20µl of the supernatant are injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30°C using a reversed phase column in two runs. One run lasts 10 minutes. The chromatograms are recorded by a fluorescence detector. The quantification is performed with the delivered EDTA-blood calibrator; the concentration is calculated via integration of the peak height by the external standard method for the reduced fraction and the internal standard method for the total GSH fraction. For the measurement of the reduced fraction it is not possible to use the internal standard because of the production of mixed sulfides.

The amount of oxidized glutathione can be calculated by subtraction of:

$$\text{Glutathione}_{\text{total}} - \text{Glutathione}_{\text{reduced}}$$

Note that the obtained difference must be divided by two as oxidized glutathione (GSSG) is composed of two GSH molecules.

Summary:

Besides many other parameters the advantage of HPLC method lies in the simultaneous handling of many analytes in a single test. The HPLC system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical routine determination in a quick and precise manner. Unlike immuno assays with up to six calibrators per test, a one-point calibration is mostly sufficient to calibrate the test system. It is possible to automate the sample application and calculation of the results so that even higher numbers of samples can be handled nearly without control.

4. MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
KC 1800LM	MOPHA	Mobile phase	1000 ml
KC 1800KA	CAL	Calibrator, lyophilized	8 vials
KC 1800IS	INT STD	Internal standard	6 ml
KC 1800T	RECSOL	Reconstitution solution	5 ml
KC 1800VL	DIL	Dilution solution	25 ml
KC 1800RL	REDSOL	Reduction solution (lyoph. 1,2 ml)	1 vial
KC 1800RB	REABUF	Reaction buffer	27 ml
KC 1800DL	DER	Derivatisation solution	12 ml
KC 1800FR	PREC	Precipitation solution	12 ml
KC 1800KO	CTRL 1 CTRL 2	Control 1 and 2, 250 µl lyophilized	2 x 3 vials

HPLC column (KC 1800RP) as well as individual components can be ordered separately from Immundiagnostik. Please ask for the price list of the individual components.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with Fluorescence-detector
- Reversed phase C₁₈-column
- Oven or water bath for heating at 60 °C

6. PREPARATION AND STORAGE OF REAGENTS

- Reconstitute **calibrator** (CAL) in **250 µl** reconstitution solution (RECSOL). The content of one vial is intended for one examination. Discard remaining unused material. The concentration of glutathione might have minor changes from lot to lot.
- Reconstitute **controls** (CTRL1, CTRL2) in **250 µl** reconstitution solution (RECSOL). Take the exact concentration from the product data sheet.
- Reconstitute **reduction solution** (REDSOL) in **1.2 ml** reconstitution solution (RECSOL). The reconstituted reduction solution is stable for 3 months when stored at **2-8 °C**.
- **The lyophilized reagents should be stored at -20 °C**. All other test reagents are stable at **2-8 °C**, up to the date of expiry stated on the label.

7. PRECAUTIONS

- For in vitro diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Precipitation solution contains acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water. Do not breathe vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on the kit label.

8. SPECIMEN COLLECTION AND PREPARATION

EDTA-blood is suited for this test system.

Glutathione is quite sensitive against oxidation. Samples should be transported at 2-8 °C.

Samples are stable for at least 2 days at 2-8 °C or 2 weeks at -20 °C. During longer storage, the content of oxidized glutathione increases.

Note: Before analysis, lyse the erythrocytes by freezing and thawing in order to release glutathione.

9. ASSAY PROCEDURE

Procedural notes

- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

For Research Use Only
For Reference Only

Sample preparation

Pipette in a 1.5 ml reaction tube (e.g. Eppendorf):

100 µl Patient sample, calibrator (CAL) or control (CTRL1, CTRL2)

+

200 µl dilution solution (DIL)

mix

Total-glutathione

50 µl diluted sample

+

100 µl internal standard (INT STD)

+

20 µl reduction solution (REDSOL)

+

100 µl derivatisation solution (DER)

Reduced glutathione

50 µl diluted sample

+

100 µl reactions buffer (REABUF)

+

100 µl derivatisation solution (DER)

Mix well, incubate for **20 minutes** at 60 °C

Add **100 µl** precipitation solution (PREC), mix well

Precipitate for **10 minutes** at 2 -8 °C and centrifuge for **10 min** at 10.000 g

Add **200 µl** supernatant to **200 µl** reaction buffer (REABUF) in auto sampler-vials, mix well

Inject **20 µl** in the HPLC-system

Chromatographic conditions

Column material:	MZ Inertsil ODS-2; 5 µm		
	MZ PerfectBond ODS-2; 5 µm		
	Bischoff Prontosil Eurobond; 5 µm		
Column dimension:	125 mm x 4 mm		
Flow rate:	0.75 - 1.0 ml/min		
Detection :	Fluorescence:	Excitation	385 nm
		Emission	515 nm
Temperature:	30 °C		
Injection volume:	20 µl		
Running time:	10 min		

It is recommended that a guard column is used to extend column life.

10. TREATMENT OF THE COLUMN

After analysis, the column should be flushed with 30 ml aqua bidest. (1 ml/min) and stored in 50% methanol in aqua bidest. (approx. 30 ml, flow 0.7 ml/min). Before use, the system should be equilibrated with ca. 30 ml eluent.

11. RESULTS

Calculation

Calculation of glutathione_{total}


$$\frac{\text{Peak height sample} \times \text{Concentration of the calibrator}}{\text{Peak height internal standard in the sample}} \times F = \text{Concentration sample}$$

$$F = \frac{\text{Peak height internal standard in the calibrator}}{\text{Peak height calibrator}}$$

Calculation of glutathione_{reduced}

$$\text{Concentration sample (nmol/l)} = \frac{\text{Peak height sample} \times \text{Concentration of the calibrator}}{\text{Peak height calibrator}}$$

Typical chromatogram

12. LIMITATIONS

Do not use serum or plasma samples. The content of glutathione in serum and plasma is lower, and it is not possible to distinguish between the oxidized and reduced glutathione form. Do not use lipemic samples.

13. QUALITY CONTROL

Expected values

(n=50)

GSH _{total} :	763 – 1191 $\mu\text{mol/l}$
GSH _{reduced} :	620 – 970 $\mu\text{mol/l}$
GSH _{reduced/total} :	81 – 93 %

It is recommended that each laboratory should establish its own normal range. Above mentioned values are only for orientation and may vary from other published data.

Controls

Control samples should be analyzed with each run. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient samples may not be valid, if within the same assay one or more values of the quality control sample are outside the acceptable limits.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay CV: GSH_{total}: 3.9 % [n = 12]
GSH_{reduced}: 3.3 % [n = 12]

Inter-Assay CV: GSH_{total}: 4.2 % [n = 12]
GSH_{reduced}: 3.3 % [n = 12]

Linearity

up to 10 mmol/l (GSH_{total} and GSH_{reduced})

Detection limit

6 µmol/l (GSH_{total} and GSH_{reduced})

15. DISPOSAL

The MOPHA (mobile phase), REDSOL (reduction solution), INT STD (internal standard), and DER (derivatisation solution) must be disposed as non-halogenated solvents. The PREC (precipitation solution) can be neutralized with NaOH to pH 7.0 and disposed as salt solution.

Important: Reaction will produce heat, be careful!

Please refer to the appropriate national guidelines.

16. TROUBLESHOOTING

Problem	Possible reason	Solution
No signal	No or defect connection to evaluation system	Check signal cord and connection
	Detector lamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Double peaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Auto sampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flow cell is dirty	Clean flow cell

17. REFERENCES

Henning S.M., J.Z. Zhang, R.W. McKee, M.E. Swendseid, R.A. Jacob. Glutathion blood levels and other oxidant defence indices in men fed diets low in vitamin C. American institute of nutrition, 1991, 1969-1975.

Rajender K.C., F.W. Lewis, M.H. Kutner, D.M. Bate, R.G.B. Roy, D. Rudman. Plasma cysteine, cystin, and glutathione in cirrhosis. Gastroenterology 1984; 87; 770-776

Meister A. Mitochondrial changes associated with glutathione deficiency. Biochim Biophys Acta 1995; 35-42.

18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes must be avoided.
- All reagents in the test package are for in-vitro diagnostic use only.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not interchange different lot numbers of any kit component within the same assay.
- Quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.